

Tutorial

1. Getting the Code

1. Download the source code from the provided repository or location.
2. Open the `Final.py` file using your preferred Python IDE, such as **Spyder**, **PyCharm**, or **VS Code**.

2. Running the Code

Step 1: Modify Input and Output Paths

Before running the script, modify the input and output file paths in the `main()` function in `Final.py`. Ensure the paths point to the correct image files.

Example:

```
IMRI = cv2.imread('path_to_MRI_image.png', cv2.IMREAD_GRAYSCALE) # Replace with your MRI image path
```

```
IPET = cv2.imread('path_to_PET_image.png', cv2.IMREAD_COLOR) # Replace with your PET image path
```

Step 2: Install Required Libraries

Ensure that Python 3.7 or later is installed on your system. Then, install the necessary libraries by running the following command in your terminal or command prompt:

```
pip install numpy opencv-python matplotlib scipy
```

Step 3: Run the Code

After modifying the input/output paths and installing the required libraries, press **Run** in your Python IDE to execute the `Final.py` script.

3. Explanation of Expected Outputs

The IEXDoG_COA algorithm outputs the following:

- **Fused Image (FColor):** A single image that combines:
 - **MRI structural details:** High spatial resolution and clear anatomical features.
 - **PET functional details:** Color-coded metabolic activity.
- **Optimal Parameters:**
 - α : Weight for the MRI low-frequency component.
 - β : Weight for the PET low-frequency component.

- **Log Information (Optional):**
 - Details of the optimization process, such as iterations and fitness values.

Medical Applications

The fused image enhances contrast and visualization, making it useful for:

- Tumor detection.
- Brain activity analysis.
- Surgical planning.